

# Decubitus Therapy Bed Hydrofloat Unit MK-1



Quiet pressure relief

**hanse**  
MEDIZINTECHNIK

experience the function, discover the special

# Hydrofloat Unit MK-1



Special care bed to reduce body contact pressure by means of the "dry flotation" method when the tissues no longer tolerate the pressure exerted by the patient's own weight.

## Some of its many advantages are:

- Increased blood supply to tissues and restored oxygen/nutrient supply to cells
- Reduced cost of nursing care
- Temperature control
- Minimized body contact pressure
- No patient turning required

## Examples of Application:

- Decubitus therapy and prevention
- Treatment of burns victims
- Immobile patient
- Long-term care
- Deep and extensive decubitus ulcers will heal within a short time while requiring much less nursing care

## Necessary Equipment:

- 1 litre disinfectant concentrate/conditioner for mattress filling

## Technical Data:

- Weight (empty): 117 kg
- Weight (with water filling): approx 317 kg
- Maximum patient weight: 200 kg
- Dimensions:  
(width x length x height) (920 x 2100 x 750) mm
- Mains voltage: 230 V, 50 Hz
- Heating: 350 W
- Power consumption: 1.7 Ampere

## Recommended Equipment:

- Monkey pole with triangular handle  
To assist patients with sitting up
- Fur cover made of anti-decubitus fleece  
Washable to 95°C and autoclavable for 30 minutes at 120°C
- Fitted sheet "anti-decubitus cover"  
Highly elastic, particularly absorbent and temperature-regulating.
- Washable to 60 °C, no fabric softeners, no spinning, no ironing