

Decubitus Therapy Mattress Easyfloat



Quiet pressure relief

hanse
MEDIZINTECHNIK

experience the function, discover the special

Easyfloat



The unique water floatation system of the “Easyfloat” decubitus therapy mattress offers maximum pressure relief for decubitus ulcers. The “Easyfloat” fits into all standard care beds and allows a 30 degree raised upper body.

Advantages:

- No patient turning required
- Decubitus therapy mattress with no internal noise
- Stepless temperature adjustment to prevent hypothermia
- Prone position and recovery position are possible

Special Features:

- Clean and germproof design
- Water sack made of highly flexible PVC, bacteriostatic and skin friendly
- Low power consumption

Technical Data:

- Weight (empty): 18 kg
- Load bearing: For patients of all weight, but depending on the care bed's load bearing capacity
- Width: 900 mm
- Length: 1950 mm
- Height: 160 mm
- Water filling: 50 - 65 litres
- Heating: 135W, 230V, 50Hz, adjustable from 24°C to 36°C
- Safety thermostat: Set at 37.5°C

Typical Applications:

- Long-term care
- Pain therapy
- After skin transplants
- Spinal problems
- For burns victims
- General wound healing and skin complaints
- Anorexia
- Emaciation

Necessary Equipment:

- 1 litre waterbed conditioner/disinfectant cleaner

Recommended Equipment:

- Fur cover made of anti-decubitus fleece (washable to 95°C, autoclavable for 30 minutes at 120°C)
- Fitted sheet “Anti-Decubitus Cover”

- Medical device -

