

Shower and Bathing Trolley Anna



Transport · Showering · Bathing

hanse
MEDIZINTECHNIK

EXPERIENCE THE FUNCTION - DISCOVER THE SPECIAL

Anna



- Easy patient transfer with slide mat.
- Narrow width to facilitate patient transfer.
- Fold-down sides for extra space.
- Can even be used in narrow areas.
- Space-saving storage.



- Closing the sides creates a generously sized lying area. The sturdy sides offer freedom of movement as well as security.
- Shower and bathing in up to 15 cm of water
- Rounded edges and smooth surfaces for easy cleaning and hygiene.



- The patient rests comfortably on a padded slide mat. A waterproof pillow is available to support the head.
- The shower trolley has the movement trendelenburg as well as reverse trendelenburg.



- Remote control handset for height adjustment to the carers's optimum working height or to the patient's individual bed height.
- Big, smooth-running castors and an adjustable steering castor.
- Central braking system

Technical Data:

- Lying area: 730 x 1990 mm
- Height adjustment: 634 - 934 mm (lying surface)
- Maximum width: 840 mm
- Maximum length: 2061 mm
- Maximum water depth: 150 mm
- Weight: 106 kg
- Maximum patient weight: 150 kg

Accessories:

- Pillow, waterproof
- Support for upper raised body
- Hand shower with valve, pressure hose, coupling
- Spare battery set

hanse
MEDIZINTECHNIK
 Dipl. Ing. P. Hettmer GmbH

- Medical device - **CE**

Zeiss-Str. 42 · D-23626 Ratekau
 Phone +49(0)4504/8182-0 · Fax -20
 info@hanse-medizintechnik.de
 www.hanse-medizintechnik.de

Subject to technical/design modifications 03/08