

Mobilisation-Rehabilitation-Wheelchair Thessa



Easy Transfer · Safe Transport · Good Mobilisation

EXPERIENCE THE FUNCTION – DISCOVER THE SPECIAL

hanse
MEDIZINTECHNIK

Thessa



Thessa E with tilt table position up to 48°

This compact and easy to manoeuvre chair combines wheelchair and stretcher in one unit.

The seat, back and leg sections are individually and continuously adjustable and can be combined as required. When the chair is set to the lying position, the armrests lower themselves automatically. In the sitting position, the armrests can be folded up individually to improve patient transfer. All cushions are removable for ease of cleaning.

Applications range from sitting in comfort to providing mobilisation and activation for patients highly dependent on care in hospitals, nursing homes, homes for the disabled, hospices and at home.



Lying position for easy transfer



Stand-up aid



Relaxed chair position

Two models are available:

Thessa M with manual mechanical adjustment

Thessa E which is driven by an electric motor (mains or battery powered) and also tilt table functions, up to 48°.

- Continuously adjustable from sitting to lying position
- Adjustment of seat/stretchers tilt
- Movement trendelenburg up to 16°, reverse trendelenburg of 17°
- Tilt table, up to 48° - Thessa E
- Electrically adjustable - Thessa E
- Stand-up aid of 17°
- Stable non-tip chassis
- Special smooth-running castors for high manoeuvrability
- Padded and individually adjustable armrests
- Stable footrest
- Removable cushions for optimum cleaning

Technical details:

Max width:	690 mm
Max length:	1.910 mm
Max height:	1.480 mm
Max weight: Thessa M:	65 kg
Thessa E:	70 kg
Max load:	150 kg

Detailed information see system description.