

# Toilet and Shower Chairs Wolfgang and Privat



**High stability and durability**

EXPERIENCE THE FUNCTION DISCOVER THE SPECIAL

**hanse**  
MEDIZINTECHNIK

# Wolfgang and Privat



The both models Wolfgang and Privat have an anatomically shaped seat with large hygiene opening. Smooth-running thanks to large, shock-absorbing castors and lockable castor for directional control. Robust and hard-wearing thanks to the stainless steel frame and impact and temperature-resistant plastic parts, which are individually removable. Non-tip stability. Wolfgang can be positioned over a toilet.



- Footrest is safe to stand on (Wolfgang)
- Positionable over a toilet; seat height 540 mm (Wolfgang)
- Individually removable and dirt-resistant plastic parts
- Impact and temperature-resistant plastic parts
- Robust stainless steel frame
- Anatomically shaped seat
- Large hygiene opening
- Individually adjustable armrests

#### Accessories:

- Bedpan holder
- Stainless steel bucket with lid
- Safety belt
- Padded insert to cover the opening
- Anti-decubitus seat pad with opening
- Back cushion
- Padded armrests, removable

#### Technical Data:

#### Wolfgang:

#### Privat:

Depth x width:	900 x 550 mm	710 x 550 mm
Height:	1080 mm	970 mm
Height of handle bar:	1010 mm	900 mm
Seat height:	540 mm	440 mm
Max. patient weight:	150 kg	150 kg
Diameter of castors:	125 mm	100 mm
Weight:	15 kg	13 kg



Seat cushion with opening



back cushion

